# SIEMENS

## MICROALBUSTIX®

### Reagent Strips for Determining Albumin and Creatinine in Urine

**INTENDED USE:** MICROALBUSTIX Reagent Strips are firm plastic strips that contain two reagent areas that test for albumin and creatinine in urine. An albumin-to-creatinine ratio can also be determined, which allows for the use of single-void specimens in testing. The ratio is given in milligrams albumin per gram or millimole creatinine (mg/g or mg/mmol). This product provides semi-quantitative results and can be used for screening samples for microalbuminuria; positive results should be confirmed with quantitative methods for albumin. The strips are for professional, *in vitro* diagnostic use (IVD). Read the insert carefully before using the product (□).

Test results may aid clinicians in the detection of patients at risk of developing kidney damage. Microalbuminuria has been reported to be an early predictor of the development of glomerular damage in the absence of overt nephropathy. Patients with diabetes and hypertension are the primary risk groups. Patients who have been exposed to nephrotoxins or who suffer from immune disorders form the secondary risk groups; microalbuminuria has also been reported as an early predictor of the development of pre-eclampsia during pregnancy.

#### SUMMARY AND EXPLANATION

The reagent test areas on MICROALBUSTIX Reagent Strips are ready to use upon removal from the bottle and the reagent strip is disposable. The strips are read visually, requiring no additional laboratory equipment. The directions must be followed exactly. The reagent strip must be stored in the bottle with the cap tightly closed to maintain reagent reactivity. To obtain optimal results, testing should be done on FRESH urine.

**SPECIMEN COLLECTION AND PREPARATION:** Test the urine specimen as soon as possible. Boric acid, at a concentration of 1.0 g/L, is the most effective preservative recommended. If testing cannot be done within two days after voiding, refrigerate the specimen immediately and let it return to room temperature before testing. Prolonged exposure of unpreserved urine to room temperature may result in microbial proliferation with resultant loss of albumin. Specimens may be stored at 0° to 8°C for one week or at -20°C for one month without significant effect on results with this test.

Any single-void urine specimen, when evaluated in conjunction with the albumin-to-creatinine ratio, can be used to discriminate between normal and abnormal levels of microalbuminuria. First-morning specimens are recommended. Urinary albumin fluctuates day-to-day; therefore, testing three urine samples over a three-to-six month period may increase the predictive value, where two positive samples are predictive of incipient nephropathy. Twenty-four hour or timed collections may also be used with this test to determine the albumin excretion rate (AER). Additional information on screening guidelines is available from the American Diabetes Association Position Statement.

It is suggested that you screen the urine for protein, using a Siemens Healthcare Diagnostics Reagent Strip such as MULTISTIK® 10 SG. Urines showing a result of 30 mg/dL (0.3 g/L) or greater for protein or that are visibly bloody should not be tested using MICROALBUSTIX Reagent Strips.

Contamination of the urine specimen with soaps, detergents, antiseptics, or skin cleansers may affect test results. The user should determine whether the use of such products is warranted.

**CAUTION: Ensure that work area and specimen containers are always free of detergents and other contaminating substances. Some substances can interfere with patient results.**

#### PROCEDURE: MUST BE FOLLOWED EXACTLY TO ACHIEVE RELIABLE TEST RESULTS.

1. Collect a FRESH urine specimen in a clean, dry container.
2. Remove one strip from bottle and replace the cap tightly.
3. Dip the test pads into the urine, making sure both pads are wetted.
4. Immediately remove the strip, dragging the edge of the strip against the rim of the urine container to remove excess urine.
5. Compare reagent areas to corresponding Colour Chart on the bottle label at the time specified. Read the albumin test 50 seconds after dipping and the creatinine test at 60 seconds.
6. Record the results you obtain, then discard the strip according to your standard laboratory procedures.
7. Refer to "Interpretation of Results" later in this insert to determine the albumin-to-creatinine ratio.

**INTERPRETATION OF RESULTS:** The results that can be obtained are found on the Colour Chart on the bottle label. If desired, the albumin result can be compared to the creatinine result to obtain an albumin-to-creatinine ratio (mg albumin/g creatinine or mg albumin/mmol creatinine). This ratio can be used to discriminate between normal and abnormal levels of microalbumin on single-void specimens. Use the following table to determine the ratio:

Albumin	Creatinine mg/dL (mmol/L)				
	10 (0.9)	50 (4.4)	100 (8.8)	200 (17.7)	300 (26.5)
10 mg/L	Recollect specimen*				
30 mg/L				Normal	
80 mg/L			Abnormal		
150 mg/L					

\*Specimen too dilute to accurately determine ratio result. Repeat test on new specimen, preferably a first-morning collection.

For example, if the albumin result is "30 mg/L" and the creatinine result is "200 mg/dL (17.7 mmol/L)", the albumin-to-creatinine ratio is "Normal." The ratio results are equivalent to the following values (mg albumin/g or mg creatinine):

Normal: < 30 mg/g (< 3.4 g/mmol)  
Abnormal: 30–300 mg/g (3.4–33.9 g/mmol) (indicates microalbuminuria)  
High Abnormal: > 300 mg/g (> 33.9 g/mmol) (indicates clinical albuminuria)

#### EXPECTED VALUES:

**Albumin** [ALB]: Albumin is normally present in urine at concentrations of less than 20 mg/L. Microalbuminuria is defined as an albumin excretion rate of 30–299 mg/24 hours.<sup>1</sup> Urinary albumin excretions can be temporarily elevated by exercise, urinary tract infections, and acute illness with fever.

**Creatinine** [CREATININE]: Creatinine is normally present in urine at concentrations of 10–300 mg/dL (0.9–26.5 mmol/L).

**Albumin-to-Creatinine Ratio** [ALB:CREATININE]: Albumin is normally present in urine at concentrations of less than 20 mg albumin/g creatinine (3.4 mg albumin/mmol creatinine). Microalbuminuria is indicated at a ratio result of 30–300 mg/g (3.4–33.9 mg/mmol) (Abnormal) and clinical albuminuria at a ratio result of > 300 mg/g (> 33.9 mg/mmol) (High Abnormal).<sup>2</sup>

**QUALITY CONTROL:** For best results, performance of reagent strips should be confirmed whenever a new bottle is first opened by testing with commercially-available negative and positive controls that include values for microalbumin and creatinine. These control specimens may also be randomly hidden in each batch of specimens tested. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

**STORAGE AND HANDLING:** Store at room temperature between 15–30°C (59–86°F). Do not use product after expiration date (□). Do not store the bottle in direct sunlight. All unused strips must remain in the original bottle. Transfer to any other container may cause reagent strips to deteriorate and become unreactive. Do not remove desiccant packet(s) from bottle. Remove each strip from the bottle immediately before it is to be used for testing. Replace cap immediately and tightly after removing reagent strip. Do not touch test areas of the reagent strip. Work areas and specimen containers should be free of detergents and other contaminating substances.

For visual use only

Dip test areas in urine completely, but briefly, to avoid dissolving out the reagents. Read test results carefully at the times specified, in a good light (such as fluorescent) and with the test area held near the appropriate Colour Chart on the bottle label. Do not read the strips in direct sunlight.

**NOTE:** You may notice a "burst" of colour on the albumin pad that fades within several seconds after dipping. This does not indicate a positive result, unless the colour remains until the 50 second reading time.

**IMPORTANT: PROTECTION AGAINST AMBIENT MOISTURE, LIGHT AND HEAT IS ESSENTIAL TO GUARD AGAINST ALTERED REAGENT REACTIVITY:** Discoloration or darkening of reagent areas may indicate deterioration. If this is evident, or if test results are questionable or inconsistent with expected findings, the following steps are recommended: (1) confirm that the product is within the expiry date shown on the label; (2) check performance against known positive and negative control materials; (3) restet with fresh product. If proper results are not obtained, consult your local laboratory representative for advice on testing technique and results.

**LIMITATIONS OF PROCEDURES:** As with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single result or method. The presence of haemoglobin or myoglobin (> 5 mg/dL or a visibly bloody urine) may cause falsely elevated results with both the albumin and creatinine tests. Contamination of the urine specimen with soaps, detergents, antiseptics, or skin cleaners, or the use of urine preservatives may affect test results. The presence of cimetidine may cause falsely elevated results with the creatinine test.

Substances that cause abnormal urine colour, such as drugs containing azo dyes, nitrofurantoin, and riboflavin, may affect the readability of the reagent areas on urinalysis reagent strips. The colour development on the reagent pad may be masked, or a colour reaction may be produced on the test pad that could be interpreted as a false positive.

**SPECIFIC PERFORMANCE CHARACTERISTICS:** Specific performance characteristics are based on clinical and analytical studies. In clinical specimens, the specificity of the reagent strips depends upon the variability of colour perception, the presence or absence of inhibitory factors (see LIMITATIONS OF PROCEDURES section); and the lighting conditions under which the product is used.

The performance characteristics of the MICROALBUSTIX test was determined at two hospital clinical laboratories with urine specimens from patients presenting for routine urinalysis. The albumin results were compared to the results obtained with a commercially-available immunoassay test; the creatinine results were compared to a kinetic enzyme creatinine assay.

The following table shows the results of the testing. Accuracy is defined as the agreement between the MICROALBUSTIX results and the comparative methods. Sensitivity is defined as the percentage of positive results obtained by the Reagent Strip relative to those obtained by the comparative methods, while specificity refers to the percentage of negative results.

Albumin Results:	Accuracy 89% (n = 810)	Sensitivity 90% (n = 402)	Specificity 88% (n = 408)
Albumin-to-Creatinine Ratio:	Accuracy 88% (n = 810)	Sensitivity 84% (n = 432)	Specificity 91% (n = 378)

**Albumin:** In contrived urine, the albumin test generally detects albumin at a concentration of 20–40 mg/L; because of the inherent variability of clinical urines, lesser concentrations may be detected under certain conditions. The test is specific for albumin and is not affected by the following proteins when tested at concentrations at least nine times greater than the excretion rate considered to be abnormal:<sup>3</sup>

IgG	Bence-Jones protein
IgA	Tamm-Horsfall glycoprotein
IgM	β <sub>2</sub> -microglobulin
Haptoglobin	α <sub>1</sub> -antitrypsin
Transferrin	retinol binding protein

**Creatinine:** This test area detects urinary creatinine in concentrations as low as 10 mg/dL (0.9 mmol/L); the absence of creatinine in a specimen cannot be determined.

**Albumin-to-Creatinine Ratio:** In contrived urines, the ratio generally detects 30–300 mg/g (3.4–33.9 mg/mmol) albumin-to-creatinine.

#### CHEMICAL PRINCIPLES OF PROCESSES:

**Albumin:** This test is based on dye binding using a high affinity sulfonaphthalene dye. At a constant pH, the development of any blue colour is due to the presence of albumin. The resulting colour ranges from pale green to aqua blue.

**Creatinine:** This test is based on the peroxidase-like activity of a copper creatinine complex that catalyzes the reaction of disopropylbenzene dihydroperoxide and 3,3'-dihydroxy-5,5'-tetramethylbenzidine. The resulting colour ranges from orange through green to blue.

**REAGENTS:** (based on dry weight at time of impregnation)

**Albumin:** 1.9% w/w bis (3,3'-diido-4,4'-dihydroxy-5,5'-dinitrophenyl)-3,4,5,6-tetrabromosulfonephthalene; 94.2% w/w buffer; 3.9% w/w nonreactive ingredients.

**Creatinine:** 2.5% w/w copper sulfate; 4.5% w/w disopropylbenzene dihydroperoxide; 2.0% w/w 3,3',5,5'-tetramethylbenzidine; 56.4% w/w buffer; 34.6% w/w nonreactive ingredients.

**AVAILABILITY:** Bottles of 25 strips (Catalog #2087).

#### BIBLIOGRAPHY:

1. Mogensen, C.E.: Microalbuminuria Predicts Clinical Proteinuria and Early Mortality in Maturity-Onset Diabetes. *N. Engl. J. Med.* 310: 356–360; 1984.
2. Mogensen, C.E. et al.: Prevention of Diabetic Renal Disease with Special Reference to Microalbuminuria. *Lancet* 346: 1080–1085; 1995.
3. Kaplan, K.M.: Microalbuminuria a Risk Factor for Vascular and Renal Complications of Hypertension. *Am. J. Med.* 92: 85–125; 1992.
4. Nissel, H. et al.: Renal Function in Gravidas with Chronic Hypertension With and Without Superimposed Preeclampsia. *Hypertens. Pregnancy* 15: 127–134; 1996.
5. Nathan, D.M.; Rosenbaum, C.; and Protasowicki, V.D.: Single-void Urine Samples Can Be Used to Estimate Quantitative Microalbuminuria. *Diabetes Care* 10: 414–418; 1987.
6. Gindberg, J.M. et al.: Use of Single Voided Urine Samples to Estimate Quantitative Proteinuria. *N. Engl. J. Med.* 309: 1543–1546; 1983.
7. Post, G.O. et al.: Diabetic Nephropathy. *Diabetes Care* 20: 524–527; 1997.
8. Cowell, C.T., Rogers, S., and Silkin, M.: First Morning Urinary Albumin Concentration is a Good Predictor of 24-Hour Urinary Albumin Excretion in Children with Type 1 Diabetes Mellitus. *Diabetes Care* 20: 97–99; 1997.
9. Kodama, H. et al.: The Range of Albumin Concentration in the Single-Void First Morning Urine of 1090 Healthy Young Children. *Diab. Research Clin. Proc.* 9: 55–58; 1990.
10. Manoli, R. et al.: Prevalence of Hypertension and Microalbuminuria in Adult Type 1 (Insulin-Dependent) Diabetic Patients Without Renal Failure in Italy—Validation of Screening Techniques to Detect Microalbuminuria. *Acta Diabetol.* 29: 156–166; 1992.
11. Pugia, M.J. et al.: Comparison of Urine Dipsticks with Quantitative Methods for Microalbuminuria. *Eur. J. Clin. Chem. Biochem.* 35(9): 693–700; 1997.

**TRADEMARKS:** Refer to the carton of the product you are using for the applicable Siemens trademarks.

**TECHNICAL ASSISTANCE:** For technical support, contact your local technical support provider or distributor. [www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

Siemens Healthcare Diagnostics, Sir William Siemens Square, Frimley, Camberley, Surrey, GU16 8QD, UK

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