

Vico

prompt action at all times



Haelvoet nv

Edition 03/2021

Haelvoet reserves the right to modify the product. Small colour deviations are possible.

OO1 Sideguards (*) • Safety has no price

The Atmosphere sideguards protect the patient over 2/3 of the stretcher.

002 **Mattress** • For a comfortable lying position

Two mattress types have been developed for the Vico Stretcher. The standard mattress consists of 7 cm of CMRH foam. The comfort mattress (*) has a thickness of 10 cm, of which 3 cm is filled with visco-elastic foam.

Push handles • Easy to fold

These practical handles help you to manoeuvre the stretcher. When not used, these push handles can be folded very easily. By doing so, all necessary treatments can be given without any hindrance. The push handles are standard at the head end. They can be optionally installed at the foot end.

004 **Mattress platform** • X-ray transparency

The mattress platform consists of 2 standard solid laminate sections that can be used for a radiological examination of the thorax. A mattress platform with 4 solid laminate sections (*) is optionally available.

005 **Cases** • Optimal care

Two accessory cases at the head and foot end make sure that the stretcher can be equipped optimally.

6 **Operating lever** • Within reach

The backrest can be easily adjusted up to 90° by means of a gas spring.

007 High/low system • Highly stable

The high/low supports guarantee an excellent stability at all heights. The stretcher can be adjusted between 46 cm and 85,5 cm (with castors ø 150 mm and no mattress). This simplifies the treatment of a patient.

008 Cover cap (*) • Hygienic

The removable rounded synthetic cover cap makes it easier to clean the carriage.

O09 High/low pedal • Easy to operate

The high/low adjustment is guaranteed by a hydraulic pump that can be adjusted by means of a pedal

5th castor (*) • Your helping hand

A trump when manoeuvring in small rooms and narrow corridors.

011 **Castors** • Superb manoeuvrability

Our extensive gamut of high-quality design castors guarantees a highly mobile and manoeuvrable stretcher.

O12 **Central brake** • Activated with one single movement

The castors are equipped with a central brake mechanism and accessible pedals. This simplifies patient treatments considerably.

013 (Reverse) Trendelenburg (*) • Highly ergonomic

When activating the operating pedal, a smooth inclination is created until the desired (Reverse) Trendelenburg position is obtained.

Bumpers • Nonstop protection

The Vico Stretcher is standard equipped with 4 bumpers that both protect the wall and the stretcher during any movement.

* option for the Vico Stretcher


Efficient in use

The Vico Stretcher has been designed to give first aid quickly and efficiently in several situations. Numerous possibilities and accessories make the work of the nursing staff easier, whereas the patient's welfare is increased.

Operations

The stretcher is very easy to operate and can be put in its desired position in no time.





(Reverse) Trendelenburg position

The (Reverse) Trendelenburg position (*) can be activated by means of the red pedal at the foot end of the stretcher. When the pedal is stepped on, the stretcher inclines fluently until the desired position is reached.



Seating position

The two-section mattress platform is equipped with a lever at the top of the back support. This lever is within easy reach, making it possible to obtain a comfortable patient position quickly. A four-section mattress platform is optionally available.








High/low position

The high/low position can be found and activated very fast. Thanks to the hydraulic high/low pedals on both sides of the stretcher, the mattress platform height can be adjusted ergonomically.



X-ray permeable mattress platform

The mattress platform of the Vico Stretcher is suited to take x-rays. Together with the x-ray cassette holder (*), it is possible to make a diagnosis quickly. Depending on your wishes, this cassette holder can be mounted on the backrest or put completely under the mattress platform. The patient doesn't have to be moved for this, so unnecessary risks are avoided.



Functional design

Close attention has been paid to efficient care during the development of the Vico Stretcher. Thanks to the well-considered choice of materials, as well as the use of smooth and easily accessible surfaces, this stretcher is ready to be used swiftly.







Mattress Maximum patient comfort

Two mattress types have been especially developed for the Vico Stretcher. Both mattresses are filled with fire-retardant foam. They are very durable and easy to clean. The Vico standard mattress has a thickness of 7 cm. The visco-elastic surface of the comfort mattress is 3 cm thick (*) and offers extra patient support.



The maintenance-free hydraulic pump and gas spring make this stretcher a durable investment.

Removable cover cap Hygienic

Thanks to the rounded synthetic cover cap (*), even the surfaces that are difficult to reach can be cleaned very easily.



- 001 The (Reverse) Trendelenburg pedal (*) guarantees an utmost ergonomic adjustment.
- 002 The operating lever of the backrest is within easy reach.
- 003 The high/low pedal is installed on both sides of the stretcher, and is within easy reach for the nursing staff.
- 004 The x-ray cassette holder (*) can be installed easily and quickly on both the backrest and the mattress platform frame.
- 005 The solid laminate mattress platform sections (*) are suited to perform radiological examinations, without having to move the patient.
- 006 The standard mattress has a thickness of 7 cm.
- 007 The removable cover cap (*) makes it easy to clean the stretcher.



Mobile and safe

To move a patient quickly and easily, the stretcher has to be controllable and manoeuvrable. The large castors of the Vico Stretcher and folding push handles make a valuable contribution to this.

Central brake


With one single movement the castors are put in a neutral position, or the directional castor can be activated. The Vico Stretcher can be optionally equipped with 4 brake pedals (2 extra at the head end).





4 braked castors

4 unbraked castorsThis position makes it possible to move the stretcher in all possible directions.



1 directional castor

This castor makes it easier to move the stretcher in a straight line.



5th castor

A **central fifth castor (*)** improves the manoeuvrability of the stretcher, especially in small and narrow rooms. Furthermore, the resilient qualities of this castor absorb floor bumps and avoid brusque impacts when crossing thresholds.



Sideguards

The **Atmosphere sideguards (*)** offer additional protection while moving the patient. The sideguards can be folded quickly and easily.


Folding IV rod

The **folding IV rod (*)** can be installed at the head and/or foot end of the Vico Stretcher. It is telescopically extendible, and can be tightened by means of a locking mechanism. The IV rod is equipped with 2 IV hooks. The folding IV rod (*) is available in 2 or 3 pieces.



Push handles or bed panels

The push handles guarantee that the stretcher can be moved ergonomically. They are standard at the head end, but can be optionally installed at the foot end as well. When not used, these push handles can be folded very easily. By doing so, all necessary treatments can be given without any hindrance.



Type 1: These **folding push handles** are standard at the head end. Two push handles at the foot end are optionally available.



Type 2: The optional Vela panel (*) is made of a metal bracket and a solid laminate insert panel. The Vela panel has a standard panel locking mechanism. The bracket is optionally available in a chromium-plated edition.



* option for Vico Stretcher

6



Accessories

The accessories we offer can be connected with the existing cases and accessory holders. Haelvoet makes sure that the accessories are at least as good as our stretchers, so that you can always have absolute certainty that your product has the best quality.


Monitor holder



Monitor holder with a metal frame and a solid laminate tablet.

Dimensions of the tablet: 650 x 300 mm

DIN rail



An additional chromium-plated DIN rail can be installed on both sides of the stretcher.



The inclination indicator shows the inclination angle of the mattress platform section, on which it has been put. This can be installed per pair on the backrest or on the mattress platform frame.

Folding IV rod



A telescopic folding IV rod with 2 hooks can be installed at the head and/or foot end, and has to be tightened by means of a locking mechanism. Choice between IV rod in 2 or 3 pieces.

Removable oxygen bottle holder

Fixed oxygen bottle holder



The oxygen bottle holder consists of a strong arm and a basket made of solid steel wire. It is finished with an epoxy coating or a chromiumplated layer (*). The oxygen bottle holder can be installed in the 4 corners of the stretcher.

Removable IV rod



The removable IV rod with 2 hooks can be installed in the 4 corners of the stretcher.



The fixed oxygen bottle holder is attached to the carriage. It is available with an epoxy coating or chromium-plated layer (*).

X-ray cassette holder



The x-ray cassette holder can be easily attached to the backrest platform section. It can also be put completely under the nattress platform.



Comfort mattress




The comfort mattress (*) is 10 cm thick, of which 3 cm consists of visco-elastic foam.

^{*} option for the Vico Stretcher

Vico Stretcher Summary

• Dimensions and performances


(high/low stretcher)



- width: 80 cm, length: 217,5 cm
- height variation of the mattress platform:
- \pm 46 cm \pm 85,5 cm (* with castors Ø 150 mm)
- safe working load: 250 kg
- weight (without accessories): 85 kg
- mattress dimensions: 201 x 70,5 cm, thickness: 7 cm

• Dimensions and performances

(stretcher with a fixed height)



- width: 80 cm, length: 217,5 cm
- fixed height of the mattress platform: 73 cm (with castors Ø 125 mm)
- safe working load: 250 kg
- weight (without accessories): 70 kg
- mattress dimensions: 201 x 70,5 cm, thickness: 7 cm



Haelvoet makes sure that its ecellent price-quality ratio is not at the expense of the environment. Our company considers durability — in all possible meanings — of paramount importance. All our products have a long lifespan and are delivered with a 10-year guarantee on construction errors. The Vico Stretcher is made in accordance with ISO 9001:2015 and ISO

The Vico Stretcher complies with the European Medical Device Regulation MDR 2017/745 and all other applicable European harmonised

 ${\bf EN~ISO~14971}$ Application of risk management to medical devices.

Options Vico Stretcher Standard Vico Stretcher • 4-section mattress platform with • 2-section mattress platform with solid laminate sections solid laminate sections • Integral S castors ø 150 mm, central brake (only with H/L stretcher) • Levina castors ø 125 mm, separate • Integral S castors ø 200 mm, with • 5th castor, only with Integral S castors • Extra set of brake pedals at the brakes (only with a fixed height) central brake (only with H/L stretcher) • Push handles at the head end Vela panel at the head and/or foot end Push handles at the foot end • Atmosphere sideguards • Operating lever High/low pedal • (Reverse) Trendelenburg pedal • DIN rail • Folding IV rod • Paper roll holder • 4 bumpers Cover cap

10

