UrinAX CL/MC/MALT

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Bibliografia

1. Official methods of analysis of AOAC International, 18th Ed. 2005; Horwitz, W. et al. AOAC International Gaithersburg, Md.
2. Handbook of Media for Clinical Microbiology 2nd Ed. 2006; Atlas, R.M. et al. CRC Taylor & Francis

Abkürzungen/Abréviations/Abbreviazioni/Abbreviations

Symbole/Symboles/Simboli/Symbols

Symbol/Symbole/Simbolo	Bedeutung/Signification/Significato
	Bestellnummer/Numéro d'article/ Numero di catalogo/Catalogue number
	Packungsinhalt/Contenu du coffret/ Contenuto della confezione/Contents of kit
	In vitro diagnosticum/Usage in vitro/ Per uso diagnostico in vitro/For in vitro diagnostic use
	Hersteller/Fabricant/Fabbricante/Manufacturer
	Lagertemperatur/Température de conservation/ Temperatura di conservazione/Storage temperature
	Verwendbar bis/Date de péremption/ Data di scadenza/Use before
	Chargenbezeichnung / Numéro du lot/ Codice del lotto/Lot number
	Gebrauchsanweisung beachten/Consulter le mode d'emploi/ Leggere le istruzioni per l'uso/Consult the instructions for use
	Einmalgebrauch/Usage unique/Per uso unico/single use
	CE-Kennzeichnung/CE marquage/ Marchio CE/CE marking

Field of application


Transport medium, bacterial count and selective propagation of bacteria and fungi in urine.

Package contents

- 10 UrinAX dip inoculum culture media
- 10 patient labels

Principle

UrinAX consists of a sealed container with 3 culture media:

The CLED culture medium is used to determine the total bacterial count in the urine. The MacConkey Agar is used for the selective propagation of Enterobacteriaceae. Malt agar (MALT) is used for the selective propagation of yeast, e.g. Candida albicans.

Composition

Basic composition in g/l

CLED Agar	MacConkey Agar
Peptone	8.0
Meat extract	3.0
Lactose	10
L-cystein	0.128
Bromthymol blue	0.02
pH	7.5
MALT Agar	
Malt extract	30
	pH 4.0

Storage

UrinAX should be stored at 15–25 °C. Temperature fluctuations and drafts should be avoided, and incubation should occur at 35–37 °C within 24 hours. Do not freeze. Do not use the product after the expiry date specified on the package.

Collection of urinary specimens

The urinary specimen must be inoculated into the UrinAX culture medium immediately after collection. Test results can be affected if the patient has been treated with antibiotics. In this case, the test must be performed no sooner than 48 hours after taking the last dose of medication.

Performing the test

1. Unscrew the UrinAX tube and remove the culture medium without touching the cultures.
2. Dip the culture medium into the urine mid-stream until the agar surfaces are completely covered. If the urine quantity is insufficient, pour the urine onto the agar surfaces. Allow surplus urine to run off the culture medium. Dab off the last drops of urine with filter paper.
3. Return the culture medium to the tube and close the lid.
4. Complete the patient labels and affix them to the tube.
5. Incubate the tube in an upright position in an incubator at 37 °C for 16–24 hours.

In the case of a negative result, the cultures can be incubated for a further 24 hours to detect slow-growing bacteria.


Transportation of the inoculated UrinAX

The inoculated culture medium can be transported at 15–25 °C. Incubation should occur within 24 hours.

Interpretation of results

After incubation, assess the bacteria growth on CLED Agar. Compare the bacterial count on the CLED Agar side with the following readout:

Bacterial count on CLED Agar/ml

Notes: Bacterial counts < 10⁴ bacteria/ml are usually not considered pathological. Bacterial counts between 10⁴ and 10⁵ are an uncertain result and should be repeated. Bacterial counts of > 10⁵ bacteria/ml indicate a urinary tract infection.

A mixture of different bacteria strains on the UrinAX is very probably due to contamination of the urinary specimen.

MacConkey Agar:

Assess the presence of colonies. With a positive lactose reaction, the colonies are red. With a negative lactose reaction, the colonies are colorless. Carry out biochemical or immunological identification of the bacteria.

MALT Agar:

Assess the presence of fungi growth. The growth of yeast is indicated by whitish colonies.

Quality control

Quality controls are carried out during the manufacture of every UrinAX lot. Growth characteristics of the medium can be tested with the following strains:

- E. coli ATCC 25922
- Staphylococcus aureus ATCC 25923
- Proteus mirabilis ATCC 12453
- Candida albicans ATCC 10231

Interpret the results after a 16 to 24-hour incubation period as follows:

E. coli ATCC 25922
Yellow growth on CLED medium and pink colonies on MacConkey.

Staph. aureus ATCC 25923
Growth on the CLED medium only. Yellow colonies.

Proteus mirabilis ATCC 12453
Colorless colonies on the CLED medium and MacConkey medium.

Candida albicans ATCC 10231
Whitish colonies on the CLED and MALT media.

Disposal

All laboratories are responsible for proper disposal of waste in accordance with the legal requirements.

References

1. Official methods of analysis of AOAC International, 18th Ed. 2005; Horwitz, W. et al. AOAC International Gaithersburg, Md.
2. Handbook of Media for Clinical Microbiology 2nd Ed. 2006; Atlas, R.M. et al. CRC Taylor & Francis