


FLEXIBLE DT 100 T


DT 100 T*

The extendable one.
This telescopic support
arm can be compactly
stored away following
an examination –
perfect where space is
at a premium.




*Illustrated arm shown with optional equipment

ATTRACTIVE



DT 100 T

Its functionality and design reflect the current demands on an intelligent support arm. It can be extended continuously from 84 cm to 132 cm, with a locking hinge moving vertically within a radius of 150 degrees, and with a cushioned horizontal pivoting radius of 360 degrees. The support arm can also be operated with the small DT 80 vacuum unit.



strässle & co



strässle & co
Medizintechnik GmbH
Primelweg 5
D-72461 Albstadt

Phone +49(0)7432 220186
Fax +49(0)7432 220189

info@straessle-co.de