

enbio 

Magic Filter

Quick Start

Ver.3 04.12.2020



Use distilled water
Używaj wody destylowanej
Verwenden Sie destilliertes Wasser
Utilizzare acqua distillata
Utilizar agua destilada
Folosiți apă distilată
Utiliser de l'eau distillée
Gebruik gedestilleerd water
Используйте дистиллированную воду

 **Hier Bestellen!**

 **PRAXISDIENST**
Medizinprodukte seit 1953

Manufacturer: Enbio Group AG, Eichengasse 3, 4702 Oensingen, Switzerland






NEW
NOWY
NEU
NUOVO
NUEVO
NOU
NOUVEAU
NIEUW
НОВЫЙ



HALF USED
ZUŻYTY W POŁOWIE
HALB GEBRAUCHT
USATO A METÀ
MEDIO USADO
JUMATATE FOLOSITĂ
MOITIÉ UTILISÉE
HALF VERBRUIKT
ПОЛОВИНА
ИСПОЛЬЗОВАНИЯ



BUY NEW
KUP NOWY
NEU KAUFEN
COMPRARE NUOVO
COMPRAR NUEVO
CUMPĂRĂ NOU
ACHETER NOUVEAU
KOOP NIEUW
КУПИТЬ НОВЫЙ



PRAXISDIENST
Medizinprodukte seit 1953

Hier Bestellen!