



BOSCH
Invented for life



















Technical Datasheet.

Universal, intuitive, fast.

vivalytic from Bosch - The universal diagnostic platform.

Technical Datasheet

Description	Value	Description	Value
 Display	7 inch 16:10, 1024 x 600 pixel touchscreen	 Weight	13 kg
 Operating air pressure range	850-1,100 hPa, corresponds to pressure range 0-1, 400 m above sea level	 Storage humidity	20-95 % (not condensing)
 Operating temperature	15-30 °C	 Operating humidity	30-80 % (not condensing)
 Storage temperature	-20-60 °C	 Electrical data	100-240 V~, 50/60 Hz, 160 VA
 Data transfer	Ethernet 10/100 MB, WLAN 2.4 GHz (802.11b/g/n); internal: Bluetooth v4.1, 2.4 GHz (low energy), USB 2.0	 Instrumental safety	IEC/EN 61010-1 IEC/EN 61010-2-010 IEC/EN 61010-2-101 Directive 98/79/EC
 Electromagnetic compatibility	IEC/EN 61326-2-6 RED 2014/53/EC FCC47 CFR 15	 Memory capacity	16 GB
 Dimensions	Length 400 mm Width 204 mm Height 388 mm	 Mean loudness	≤ 55 dB(A) in operating mode. Short term loudness can exceed mean loudness
 Distance/Space to the wall	~20 cm	 Socket	Use multiple sockets for countries X, Y, Z


Valytic Workflow: Good results can be this easy.




01 Scan the sample code.



02 Scan the cartridge code.



03 Fill in the sample and close the cartridge cover.



04 Insert the cartridge. Test starts automatically.

 **Order here!**

 **PRAXISDIENST**
Medical Supplies since 1953

Bosch Healthcare Solutions GmbH Stuttgarter Str. 130, 71332 Waiblingen, Germany

www.bosch-valytic.com

Subject to technical modifications.

Printed in Germany. Any colour deviations are due to restrictions in the printing process. Printed on 100 % chlorine-free bleached cellulose, produced in an environmentally friendly process. Not all products are available in all regions. Ask your local sales representative for availability in specific markets. For in vitro diagnostic use only.