



Do for your think



# UT100V

## Veterinary Handheld Pulse Oximeter



UT100V veterinary pulse oximeter is a compact, light easy for carrying and handling device, measuring the animal's SpO2 and pulse rate. Mainly used in vet clinics and hospitals, suitable for dog, cat and horse.

[Order here!](#)

**PRAXISDIENST**  
Medical Supplies since 1953



Do for your think

### UTECH CO., LTD.

3rd floor, No. 39, Jinqiao Road, Zhangjiawan, Shapingba District, Chongqing, China  
Tel: 86-23-61728390 Fax: 86-23-61728391 www.chinautech.com lily@chinautech.com

# UT100V

## Veterinary Handheld Pulse Oximeter

### Product Highlight

- Easy to adjust LCD display, brightness to optimize viewing.
- Two selecting working modes: spot-check and monitor.
- Visual and audible alarm signals.
- Trend data storage time up to 120 hours.
- Trend data can be transferred to PC for printing, analysis and viewing.
- Functional of management data analysis.
- Veterinary special SpO2 sensor.



### Technical Specifications

#### SpO2

range..... 0-100%  
 accuracy..... ±2 at 70-100%  
 resolution..... 1%

#### Pulse rate

range..... 30-400bpm  
 accuracy..... ±2 at 30-400bpm  
 resolution..... 1bpm

#### Electrical

AC Power Requirements..... 100-240VAC, 47-63Hz  
 Power Consumption..... 20VA

#### Batteries

Type..... Alkaline/Ni-Mh Batteries  
 Operating Time..... >20Hours  
 Charge Time..... <6Hours

#### Environment

Operating Temperature..... 0-45°C  
 Storage Temperature..... -20-60°C  
 Operating Humidity..... 30-95%  
 Operating Altitude..... -500-5000m

### Physical Characteristics

Dimensions..... (13.5cm×7.5cm×2.8cm)

Weight..... 258g with batteries

#### Alarms

Audible and visual alarms for high and low saturation(0-100%), pulse rate(30-250bpm), Sensor off, full storage and low battery

Alarm Volume..... 85dB



### Product Accessories



SpO2 Sensor



CD



AC Charger



Alkaline batteries



Hangingstrip



## UTECH CO., LTD

3rd floor, No. 39, Jinqiao Road, Zhangjiawan, Shapingba District, Chongqing, China  
 Tel: 86-23-61728390 Fax: 86-23-61728391 www.chinautech.com lily@chinautech.com