


EN

BRYGGA®

MAVIG




Undereath the Apron

Scan QR code  
to view the web video



MAVIG website:  
[www.mavig.com/product/ra-brygga/](http://www.mavig.com/product/ra-brygga/)



**Model:** BRYGGA

Carbon brackets, weight relief for the shoulder when wearing X-ray protective aprons  
Structure: 5-piece – Material: Robust, 100% carbon fibre  
Foam padding: Exterior: Evazote – Interior: Polyethylene

Made in Sweden


## BRYGGA® – Shoulder Brace for Individual Weight Relief

Wearing a BRYGGA shoulder brace under the protective X-ray apron is an easy-to-use and efficient way to reduce the weight load on the shoulders and prevent muscle fatigue. In this unique way, the patented invention from Sweden reduces the pressure on the shoulder area that is normally caused by the weight of the apron worn.

**Optimum weight distribution:** BRYGGA distributes the weight more evenly and relies primarily on the bone structure rather than the muscles and nerve cords. This proves to be particularly advantageous for long periods of wear.

**Modular structure:** The multi-part construction allows for a quick adjustment to the respective user and ensures a pleasant and comfortable wearing feeling in an easy and uncomplicated way.

- Patented solution – Made in Sweden
- For MAVIG apron models RA631 and RA632
- Two sizes: S/M and L, individually adjustable
- No impairment for the user, ideal for long wearing times
- Promotes an ergonomically correct posture
- Significantly reduces the risk of tension and numbness in the arms
- Full freedom of movement for neck and arms
- Easy cleaning



**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)